



**Betriebsanleitung
Instruction manual
Mode d'emploi
Instrucciones de Empleo
Istruzioni per l'uso**

ROTACOOOL

CE

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EMV-Richtlinie:
EN 50051-2; EN 50082-2

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E

Thank you for buying this appliance. You have acquired a product that was manufactured and tested by Heidolph in accordance with the highest quality standards. We hope it will enable you to carry out your work easily and to your complete satisfaction.

CONTENTS

Illustrations	2
Packing list	12
Outline description	13
General information	13
Safety instructions	13
Assembly	14
Operation	16
Routine cleaning and maintenance	16
Dismantling, transport and storage	17
Disposal	17
Troubleshooting	17
Technical data	18
Warranty, exclusion of liability, copyright	18
Questions / repair work	19
Declaration of CE conformity	20

Important information

Advice about power cord / mains supply

Caution: mandatory action

PACKING LIST

Please make sure the package contains the following parts:

Description	Quantity	Order Number
Circulating chiller ROTACOOOL 230V 50Hz	1	591-00200-00
or		
Circulating chiller ROTACOOOL 115V 60Hz	1	591-00200-01
Connection for the cooling liquid M 16 x 1/NW 7 (A)	2	11-300-005-23
Cover I	1	23-30-01-03-12
Connection cable for Europe (230V 50Hz)	1	14-007-003-81
or		
Connection cable for USA (115V 60Hz)	1	14-007-003-89
Operating instructions	1	01-005-002-77

OUTLINE DESCRIPTION

L-shaped, CFC-free circulating chiller for HEIDOLPH's rotary evaporator LABOROTA 4000 ff. To save space, the rotary evaporator (P) and HEIDOLPH's vacuum pump ROTAVAC vario control (Q) can be set up on the circulating chiller. The circulating chiller is equipped with an air-cooled refrigerating machine and a circulating pump. The heat exchanger (condenser), pump and housing are made of stainless steel. The cooling liquid is contained in an atmospherically open expansion vessel while the external coolant circuit is closed. Temperature controller with digital display for actual and target values, constant control, temperature stability 0.5 K.

GENERAL INFORMATION

Unpack the appliance and its components carefully and compare the contents with the packing list. If you discover any damage or deviations from the packing list, notify your supplier immediately.

Read the instruction manual carefully and make sure that anyone else who will be operating the appliance also reads it thoroughly before use.

Keep the instruction manual accessible to everyone at all times.

This appliance is equipped as standard with a Euro plug (DIN 49441 CEE 7/VII 10/16 A 250 V) or with a US standard plug (NEMA Pub. No. WDI-1961 ASA C 73.1 – 1961 page 8 15 A 125 V) for use in North America. Where connection to a different mains supply is required you must use an approved adapter or arrange for an electrician to replace the supplied plug with a 3-pin plug designed and approved for the mains supply in question. The appliance is delivered in earthed condition. If the original plug is ever changed, make sure that the protective earthed conductor from the appliance is connected to the new plug. Color code used for wiring inside the appliance:

GREEN/YELLOW = earth (PE)
 BLUE = neutral (N)
 BROWN = live (L)

IN GREAT BRITAIN, AN APPROVED STANDARD 13A PLUG TO BS 1363 (3-PIN PLUG) IS NECESSARY.

SAFETY INSTRUCTIONS

Follow the safety regulations in force for your laboratory!

Be extremely careful when dealing with inflammable substances.

Before you connect the appliance to your power supply, make sure that the supply voltage is the same as that on the rating plate.

E

Make sure the appliance stands securely!

The appliance is not explosion-proof. Caution must therefore be exercised when using it near flammable and explosive substances.

ASSEMBLY

1. Important

Read the following instructions before switching on. Have any work required on the electrical, electronic and refrigerating components performed only by suitably qualified professionals. Disconnect from the mains before servicing!

2. Ambient conditions

Set up the appliance on a flat, horizontal surface. Make sure there are no light particles (polystyrene, paper or the like) in the area of the air intake in the base of the housing (U). The admissible ambient temperature is between + 15°C and 32°C. Protect the appliance against moisture.

The rotary condenser has to be set up with a clearance of at least 8 cm from any rear wall in order not to obstruct air escaping at the back of the unit.

3. Operating temperature range

ROTACOOOL – 10°C to + 40°C

4. Electrical connection

Switch off the ON/OFF switch (H). Observing the technical data on the type plate, connect the appliance to the mains voltage from the appliance socket (S) using the connection cable supplied.

5. Cooling liquid

Use water, ethanol or a suitable thermofluid as cooling liquid. Make sure the thermofluid's maximum and minimum admissible temperatures are suitable for the operating temperature range desired.

Use cooling liquids with a freezing point at least 5 K below the operating temperature. If you want to set an operating temperature of e.g. –10°C, the cooling liquid must be guaranteed not to freeze at temperatures down to –15°C. If water is used as cooling liquid at temperatures below +5°C, you must either add a standard anti-freeze agent or use a different cooling liquid (e.g. Ethanol) instead.

Observe the following information for thermofluids:

Maximum operating temperature: The flash point of the thermofluid used should be at least 5 K higher.

Minimum operating temperature: The viscosity of the thermofluid increases at low temperatures. This impairs flow and reduces heat transfer. Viscosity should therefore be as low as possible and should not exceed 120 mm²/s.

Remember that many thermofluids are hygroscopic. The air humidity to collect in the thermofluid will form a layer of ice on the surface of the evaporator which will impair heat transfer. At temperatures below room temperature in particular, the filling apertures should be sealed with the cover I. The water to accumulate in the thermofluid must be evaporated off from time to time, if necessary.

The following substances may not be used in the bath:

- **Cooling liquids which contain additives such as ethers, esters and amines (which are contained in some glycolenes)**
- **Demineralized or distilled water**
- **Mineral or sea water**
- **No CaCl₂ brines as these destroy the high-grade steel material 1.4301.**

The bath for the cooling liquid for the circulating chiller has a capacity of approx. 2.5 l.

6. Tube connections


Screw the connections for the cooling liquid (A) to the threaded sockets (B) using an open-ended, 19 mm spanner. To connect the condenser of the rotary evaporator to your circulating chiller, you will need a tube with a nominal width of 7 – 8 mm for both the advance and the return. Connect the tube for the advance to the left socket (D) on top of the appliance and make sure it is leak-proof. Connect the return to the right socket (E). Protect both tubes against inadvertent slippage using hose clamps. If you want to work with advance temperatures in excess of 5 to 10 K either above or below the ambient temperature, these connecting tubes must be thermally insulated as well.

To ensure the proper circulation of the cooling liquid, make sure there are no kinks in the connecting tubes.

7. Filling the cooling liquid


Only an externally closed system can be connected to the rotary chiller. Shut off the discharge valve (F) (so that the slit is vertical). After removing the cover I, fill the system through the filling aperture up to approx. 1 cm below the overflow hole (T) in the filling gauge (G). Aerate the condenser of the rotary evaporator by switching on the ROTACOOOL at the On/OFF switch (H) (press the switch and the control lamp will light up to confirm operation). The circulating condenser aerates itself via the built-in expansion vessel.

During operation of your circulating chiller and when using a large capacity condenser (e.g. condenser G4) the level of liquid in the filling gauge (G) may fall so far that the cooling liquid level switch responds and stops the circulating chiller. Top up with the cooling liquid. (See Chapter 8 „Discharging the cooling liquid“).

8. Discharging the cooling liquid

Remove the screw-on cap (I) of the discharge socket (J) and slide a tube with an inner diameter of 8 mm over the discharge socket (J). Secure the tube to prevent inadvertent slippage. Hold the other end in a collection vessel and open the discharge valve (F) (so that the slit is horizontal).

Every time before disconnecting the condenser from the rotary evaporator you must first empty the circulating chiller in the way described if a large-capacity condenser (e.g. condenser G4) is being used. This is necessary because the circulating volume is greater in this case than the capacity of the expansion vessel.

9. Setting up the rotary evaporator

The rotary evaporator (P) is set up on the flat surface of the circulating chiller. The control panel and right edge of the rotary evaporator LABOROTA 4000 ff must be flush with the edges of the circulating chiller.

10. Setting up the vacuum pump ROTAVAC vario control

The vacuum pump ROTAVAC vario control (Q) is set up on the circulating chiller. The rubber feet of ROTAVAC vario control is placed in special cups (V) on the top of the circulating chiller to prevent inadvertent slippage.

OPERATION

Before commencing operation, make sure all the tubes have been properly connected and are leak-proof as the circulating pump will start up as soon as the ON/OFF switch is actuated.

Press the ON/OFF switch (H) to commence operation. The integrated control lamp will confirm operation. Tempering and the circulating pump will start up. When the appliance is switched on, all five segments of the temperature display (K) will light up for approx. 5 seconds. After this, the actual temperature will be displayed, which means that the controller is now ready for operation. The red LED (L) will show the actual status of the refrigerating machine:

LED on: performance adjustment on
LED off: performance adjustment off

The ROTACOOOL is equipped with a level switch which means that the appliance will not start up if there is not enough cooling liquid or no cooling liquid at all in the bath (temperature display (K) not active). If this is the case, top up the cooling liquid (see ASSEMBLY).

Target advance temperature

Call up the target temperature by depressing the „SET“ key (M) next to the display. Keeping this key depressed, you can adjust the target temperature by pressing the cursor keys “▲“ (N) (which increases the value) and “▼“ (O) (which lowers the value). You will see the advance temperature again as soon as you release the “SET“ key.

ROUTINE CLEANING AND SERVICING

The cooling air for the refrigerating machine is sucked in through the liquefier in the base of the housing. The liquefier should be cleaned from time to time. To do this, tip the discharged appliance onto its back. The fins of the black liquefier (U) can now be cleaned using a brush

or vacuum cleaner. Never use pointed or sharp-edged tools for this purpose, however, make sure the fins of the liquefier are not damaged or deformed as this will impair air flow.

Use a damp cloth (mild soap solution) to wipe clean the housing and surfaces of the appliance.

IMPORTANT: Never use chlorine bleach or chlorine-based detergents, scouring agents, ammonia, cleaning wool or cleaning agents with metallic components to clean the appliance, as these will damage the surfaces of the appliance.

The appliance is maintenance-free. Any repairs which may be necessary must be carried out by an authorized Heidolph technician. Please contact your local HEIDOLPH dealer or representative to this end.

DISMANTLING, TRANSPORT AND STORAGE

Dismantling

Switch off the appliance and disconnect from the mains. Remove the rotary evaporator and vacuum pump.

Drain off all the cooling liquid still in the appliance (see ASSEMBLY).

Transport and storage

Pack the appliance and parts in the original packing and make sure they are packed upright, as indicated by the arrows on the packing. If the original packing is no longer available, make sure adequate alternative is found and that the appliance is transported upright (mark the packing with arrows).

The ROTACOOOL may not be tipped or laid on its side during transport.

DISPOSAL

Please make sure that any old appliances or defective parts are disposed of only by specialist builders of refrigeration plant (environmental disposal of the refrigerant). Please also dispose of the packing materials in an environment-friendly manner (sorting).

TROUBLESHOOTING

Have any work required on the electrical, electronic and refrigerating components performed only by suitably qualified professionals.

If the appliance does not function (the temperature display does not light up):

- Check the mains cable
- Check the ON/OFF switch

E

- The level switch has responded. There is not enough cooling liquid; top up the cooling liquid up to approx. 1 cm below the overflow hole (T) in the filling gauge.

Choice of error codes displayed:

- F1H: sensor broken, sensor resistance over limit
- F1L: sensor short circuit, sensor resistance below limit
- F0-: wrong actual value
- F5-: wrong nominal entry in memory
- EP-: EEPROM error in memory

If one of these error codes appears, give Heidolph Instruments a call (+49-9122-992068) or get in touch with your local Heidolph Instruments dealer

TECHNICAL DATA

ROTACOOOL Circulating chiller

Operating temperature range	[°C]	- 10 to + 40
Resolution/accuracy	[K]	0,1 / 0,5
Temperature stability	[K]	+/- 0,5
Cooling performance in W at an ambient temperature of 20°C		at 15°C 400 W at 0°C 300 W at -10°C 270 W
Cooling liquid capacity	[l]	2,5
Refrigerant in the chiller (provided by manufacturer)		R 134 a (CFC free) 0,2 kg
Pump's maximum delivery rate	[l/min]	10
Pump's maximum delivery pressure	[bar]	0,3
Connections for the cooling liquid	[mm]	M 16 x 1 / NW 7
Material of the heat exchanger, pump and housing		stainless steel 1.4301
Basic area (W x D)	[mm]	470 x 580
Overall height/height of flat surface	[mm]	405 / 95
Flat surface for rotary evaporator (W x D)	[mm]	470 x 405
Mains voltage		230V/50Hz or 115V/60Hz
Power rating	[W]	600
Weight	[kg]	35
Protection class		IP 20
Admissible ambient conditions		15-32°C bei 80% rel. Humidity

WARRANTY / EXCLUSION OF LIABILITY / COPYRIGHT

Warranty

Heidolph Instruments guarantees warranty that the present product shall be free from defects in material (except wear parts) and workmanship for 3 years from the date shipped off the manufacturer's warehouse.

Transit damage is excluded from this warranty.

To obtain such warranty service, contact Heidolph Instruments (Tel.: +49 - 9122 - 9920-68) or your local Heidolph Instruments Dealer. If defects in material or workmanship are found, your item will be repaired or replaced at no charge.

Misuse, abuse, neglect or improper installation are not covered by this warranty guarantee.. Alterations to the present warranty guarantee need Heidolph Instruments' consent in writing.

Exclusion Clause

Heidolph Instruments cannot be held liable for damage from improper use or misuse. Remedy for consequential damage is excluded.

Copyright

Copyright in pictures and wording of the present Instruction Manual is held by Heidolph Instruments.

QUESTIONS / REPAIR WORK

If any **aspect** of installation, operation or maintenance remains unanswered in the present Manual, please contact the following address:

For repair services please call Heidolph Instruments (phone: +49 - 9122 - 9920-68) or your local, authorized Heidolph Instruments Dealer.

Note

You will receive approval for sending your defective item to the following address:

Heidolph Instruments GmbH & Co. KG
Lab Equipment Sales
Walpersdorfer Str. 12
D-91126 Schwabach / Germany
Tel.: +49 – 9122 - 9920-68
Fax: +49 – 9122 - 9920-65
E-Mail: Sales@Heidolph.de

Note

If you are based in the United States of America, please contact Heidolph US:

Heidolph Instruments LLC
Lab Equipment Sales
2615 River Rd.
Cinnaminson, NJ 08077
Phone: 856-829-6160
Fax: 856-829-7639
E-Mail: heidolph@snip.net

Safety Information

When shipping items for repair that may have been contaminated by hazardous substances, please:

- advise exact *substance*
- take proper *protective* measure to ensure the safety of our receiving and service personnel
- *mark* the pack IAW Hazardous Materials Act

CE DECLARATION OF CONFORMITY

We declare that this product conforms to the following standards or published documents:

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© Heidolph Instruments GmbH & Co. KG